


STRAIGHT AESTHETIC ABUTMENTS

Standard platform

ANGLE: 0° DIAMETER:5.0



DIAMETER	K	GINGIVAL HEIGHT	PRODUCT CODE
5 mm	7 mm	1	DXL-10010005-S
5 mm	7 mm	2	DXL-1001010-S
5 mm	7 mm	3	DXL-1001030-S

Straight Aesthetic titanium abutment for standard platform 2.24 mm with variety dimensions gingival height. Prosthetic screw in each set. It is installed using universal or hand prosthetic insertion driver.


MATERIAL Titanium Grade 5

INTERNAL HEX
2.24

AESTHETIC ANGLED ABUTMENTS

Standard platform

ANGLE: 7° DIAMETER:5.0



DIAMETER	K	GINGIVAL HEIGHT	PRODUCT CODE	ANGLE
5.0 mm	7 mm	0,5	DXL-10010005	7°
5.0 mm	7 mm	1	DXL-10010010	7°
5.0 mm	7 mm	2	DXL-10010020	7°
5.0 mm	7 mm	3	DXL-10010030	7°
5.0 mm	7 mm	4	DXL-10010040	7°

Aesthetic Angled titanium abutment for standard platform 2.24 mm with variety dimensions gingival height. Prosthetic screw in each set. It is installed using universal or hand prosthetic insertion driver.

MATERIAL Titanium Grade 5

SCREW DRIVER
1.4HEX

Each abutment set includes




PROSTHETIC SCREW

AESTHETIC ANGLED ABUTMENTS

Standard platform

DXL
MEDICAL

ANGLE: 15° DIAMETER:5.0



DIAMETER	K	GINGIVAL HEIGHT	PRODUCT CODE	ANGLE
5.0 mm	7 mm	0,5	DXL-10020005	15°
5.0 mm	7 mm	1,5	DXL-10020015	15°
5.0 mm	7 mm	2,5	DXL-10020025	15°
5.0 mm	7 mm	3,5	DXL-10020035	15°
5.0 mm	7 mm	5	DXL-10020045	15°


MATERIAL Titanium Grade5

INTERNAL HEX
2.24

SCREW DRIVER
1.4HEX

AESTHETIC ANGLED ABUTMENTS

ANGLE: 25° DIAMETER:5.0



DIAMETER	K	GINGIVAL HEIGHT	PRODUCT CODE	ANGLE
5.0 mm	7 mm	0,5	DXL-10040005	25°
5.0 mm	7 mm	1,5	DXL-10040015	25°
5.0 mm	7 mm	2,5	DXL-10040025	25°
5.0 mm	7 mm	3,5	DXL-10040035	25°
5.0 mm	7 mm	5	DXL-10040045	25°

MATERIAL Titanium Grade5


Each abutment
set includes



Aesthetic Angled titanium abutment for standard platform 2.24 mm with variety dimensions gingival height.
Prosthetic screw in each set. It is installed using universal or hand prosthetic insertion driver.

AESTHETIC ANGLED ABUTMENTS

Narrow Platform




ANGLE: 7° DIAMETER:4.0



DIAMETER	GINGIVAL HEIGHT	PRODUCT CODE	ANGLE
4.0 mm	0,5	DXL-10010005-N	7°
4.0 mm	1	DXL-10010010-N	7°
4.0 mm	2	DXL-10010020-N	7°
4.0 mm	3	DXL-10010030-N	7°

 INTERNAL HEX
2.1

MATERIAL Titanium Grade 5




ANGLE: 15° DIAMETER:4.0



DIAMETER	GINGIVAL HEIGHT	PRODUCT CODE	ANGLE
4.0 mm	0,5	DXL-10020005-N	15°
4.0 mm	1	DXL-10020010-N	15°
4.0 mm	2	DXL-10020020-N	15°
4.0 mm	3	DXL-10020030-N	15°

SCREW DRIVER 1.2HEX

 Eachabutment
 set includes

 PROSTHETIC
 SCREW


ANGLE: 25° DIAMETER:4.0



DIAMETER	GINGIVAL HEIGHT	PRODUCT CODE	ANGLE
4.0 mm	0,5	DXL-10040005-N	25°
4.0 mm	1	DXL-10040010-N	25°
4.0 mm	2	DXL-10040020-N	25°

Angulated anatomic titanium abutment for narrow platform 2..1 mm with variety dimensions gingival height or without it. Prosthetic screw in each set. It is installed using universal or hand prosthetic insertion driver.

Attention: Ø3.3 and Ø3.5 mm diameter DXL implants have the narrow platform.