


STRAIGHT ABUTMENTS

Standard platform

DXL
MEDICAL

Ø 5.0
SHORT




DIAMETER	K	GINGIVAL HEIGHT	PRODUCT CODE
5.0 mm	4 mm	1	DXL-7010501
5.0 mm	4 mm	2	DXL-7010502
5.0 mm	4 mm	3	DXL-7010503
5.0 mm	4 mm	4	DXL-7010504
5.0 mm	4 mm	5	DXL-7010505
5.0 mm	4 mm	6	DXL-7010506

MATERIAL Titanium Grade5

STRAIGHT ABUTMENTS

Ø 6.5
SHORT



DIAMETER	K	GINGIVAL HEIGHT	PRODUCT CODE
6.5 mm	4 mm	1	DXL-7016501
6.5 mm	4 mm	2	DXL-7016502
6.5 mm	4 mm	3	DXL-7016503

MATERIAL Titanium Grade5

Titanium straight abutment for standard platform 2.24 mm .
Prosthetic screw in each set. It is installed using universal or hand prosthetic insertion driver.

Each abutment set includes




PROSTHETIC
SCREW

**SCREW DRIVER
1.4HEX**

STRAIGHT ABUTMENTS


Standard platform



DIAMETER	K	GINGIVAL HEIGHT	PRODUCT CODE
5.0 mm	7 mm	1	DXL-7020501
5.0 mm	7 mm	2	DXL-7020502
5.0 mm	7 mm	3	DXL-7020503
5.0 mm	7 mm	4	DXL-7020504
5.0 mm	7 mm	5	DXL-7020505
5.0 mm	7 mm	6	DXL-7020506

MATERIAL Titanium Grade5

STRAIGHT ABUTMENTS



DIAMETER	K	GINGIVAL HEIGHT	PRODUCT CODE
6.5 mm	7 mm	1	DXL-7026501
6.5 mm	7 mm	2	DXL-7026502
6.5 mm	7 mm	3	DXL-7026503

MATERIAL Titanium Grade5

INTERNAL HEX

2.24

Each abutment set includes




PROSTHETIC SCREW

SCREW DRIVER 1.4HEX

STRAIGHT ABUTMENTS

Standard platform

DXL
MEDICAL



DIAMETER	K	GINGIVAL HEIGHT	PRODUCT CODE
3.7 mm	5 mm	0	DXL-7010500
3.7 mm	7 mm	0	DXL-7020500
3.7 mm	9 mm	0	DXL-7030500

Each abutment set includes




PROSTHETIC SCREW

INTERNAL HEX
2.24

SCREW DRIVER 1.4HEX

MATERIAL Titanium Grade 5

Narrow platform



DIAMETER	K	GINGIVAL HEIGHT	PRODUCT CODE
4 mm	5.5 mm	1	DXL-7020401-N
4 mm	5.5 mm	2	DXL-7020402-N
4 mm	5.5 mm	3	DXL-7020403-N

Each abutment set includes



PROSTHETIC SCREW

INTERNAL HEX
2.1

SCREW DRIVER 1.2HEX


MATERIAL Titanium Grade 5

Titanium straight abutment for narrow platform 2.1 mm .
Prosthetic screw in each set. It is installed using universal or hand prosthetic insertion driver.

STRAIGHT AESTHETIC ABUTMENTS

Standard platform

ANGLE: 0° DIAMETER:5.0



DIAMETER	K	GINGIVAL HEIGHT	PRODUCT CODE
5 mm	7 mm	1	DXL-10010005-S
5 mm	7 mm	2	DXL-1001010-S
5 mm	7 mm	3	DXL-1001030-S

Straight Aesthetic titanium abutment for standard platform 2.24 mm with variety dimensions gingival height. Prosthetic screw in each set. It is installed using universal or hand prosthetic insertion driver.


MATERIAL Titanium Grade 5

INTERNAL HEX
2.24

AESTHETIC ANGLED ABUTMENTS

Standard platform

ANGLE: 7° DIAMETER:5.0



DIAMETER	K	GINGIVAL HEIGHT	PRODUCT CODE	ANGLE
5.0 mm	7 mm	0,5	DXL-10010005	7°
5.0 mm	7 mm	1	DXL-10010010	7°
5.0 mm	7 mm	2	DXL-10010020	7°
5.0 mm	7 mm	3	DXL-10010030	7°
5.0 mm	7 mm	4	DXL-10010040	7°

MATERIAL Titanium Grade 5

Aesthetic Angled titanium abutment for standard platform 2.24 mm with variety dimensions gingival height. Prosthetic screw in each set. It is installed using universal or hand prosthetic insertion driver.

SCREW DRIVER
1.4HEX

Each abutment set includes



PROSTHETIC SCREW