



HEALING CAPS


Standard Platform




Ø 4.7




HEALINGCAPS




Ø 4.8



HEALINGCAPS



Ø 6.0



DXL Medical 22

PRODUCT CODE

DXL-4444720

DXL-4444730

DXL-4444740

DXL-4444750

DIAMETER

4.7 mm

4.7 mm

4.7 mm

4.7 mm

GINGIVAL HEIGHT

3.5 mm

4.5 mm

5.5 mm

6.5 mm

MATERIAL Titanium Grade 5

SCREW DRIVER
1.4 HEX

PRODUCT CODE

DXL-4004510

DXL-4004520

DXL-4004530

DXL-4004540

DXL-4004550

DXL-4004560

DIAMETER

4.8 mm

4.8 mm

4.8 mm

4.8 mm

4.8 mm

4.8 mm

GINGIVAL HEIGHT

1 mm

2 mm

3 mm

4 mm

5 mm

6 mm

MATERIAL Titanium Grade 5

SCREW DRIVER
1.4 HEX

PRODUCT CODE

DXL-4006020

DXL-4006030

DXL-4006040

DXL-4006050

DXL-4006060

DIAMETER

6.0 mm

6.0 mm

6.0 mm

6.0 mm

6.0 mm

GINGIVAL HEIGHT

2 mm

3 mm

4 mm

5 mm

6 mm

MATERIAL Titanium Grade 5

SCREW DRIVER
1.4 HEX


HEALING CAPS

Standard platform



PRODUCTCODE	DXL-4008030	DXL-4008050
DIAMETER	8.0 mm	8.0 mm
GINGIVAL HEIGHT	3.0 mm	5.0 mm

Ø 8.0



MATERIAL Titanium Grade 5


Narrow platform

Attention: Ø 3.3 and Ø 3.5 mm diameter DXL implants have the narrow platform.

Ø 4.1



PRODUCTCODE	DXL-4004120-N	DXL-4004130-N	DXL-4004150-N
DIAMETER	4.1 mm	4.1 mm	4.7 mm
GINGIVAL HEIGHT	3.5 mm	4.5 mm	5.5 mm




MATERIAL Titanium Grade 5

COVER SCREW

Standard Platform

SCREW DRIVER
1.4HEX



PRODUCTCODE

DXL-3000001 DXL-3000021 DXL-3000100

DIAMETER

3.1mm 3.1mm 6.0mm

GINGIVAL HEIGHT


0 mm 2mm 1 mm

MATERIAL Titanium Grade5

It is possible to order cover screw as sterile. Add -S to the relevant article in your order.

Narrow Platform

SCREW DRIVER
1.2HEX



Ø 2.6



PRODUCTCODE

DXL-3000001-N DXL-3000021-N

DIAMETER

2.6 mm 2.6mm

GINGIVAL HEIGHT

0 mm 2 mm

MATERIAL Titanium Grade 5

Attention: Ø 3.3 and Ø 3.5 mm diameter
 DXL implants have the narrowplatform.

It is possible to order cover screw as sterile. Add -S to the relevant article in your order.